

# SEF<sup>2</sup> 2-Part Mammography Developer

***Optimal  
Performance  
for All Single  
and Double  
Emulsion Films***

More detail and more diagnostic information – that's what makes **SEF<sup>2</sup>** the best choice for your radiology department. In addition, to increased detail, **SEF<sup>2</sup>** can help reduce radiation exposure to patients by up to 25%.

**SEF<sup>2</sup>** developer was designed with mammography patients in mind, utilizing a single formulation that allows the development of all major film brands.

In Rapid Access applications, developer temperature may be lowered to 97 degrees from 101 degrees, helping maintain QC even further.



CBM20 Chemical  
Mixing Station



## SEF<sup>2</sup> 10 Gallon Carton

Stock Number

S60-SEF2

*2 bottles, 5 quart part A developer*

*2 bottles, 16 ounce part B developer*



www.hrsimon.com

**H.R. Simon and Company, Inc.**  
3515 Marmenco Court  
Baltimore, Maryland 21230 • USA

TOLL FREE 800-638-9460  
BALTIMORE 410-636-5555  
FAX 410-636-4415  
E-MAIL sales@hrsimon.com

S60-SEF2		SEF <sup>2</sup> Developer 10 Gallon Pack			31 lbs / 14.1 kg	
TO MAKE		USE THIS AMOUNT OF WATER			USE THESE UNITS	
Gallons	Liters	Gallons +	Ounces	Liters	Part A	Part B
5	19	3	80	13.7	1	1
10	38	7	32	27.5	2	2
S60-SEF2D		SEF <sup>2</sup> Developer 200 Gallon Bulk Mix			609 lbs / 277kg	
TO MAKE		USE THIS AMOUNT OF WATER			USE THESE UNITS	
Gallons	Liters	Gallons +	Ounces	Liters	Part A	Part B
200	757	145	0	548.8	1	1

\*Never leave containers half-full

PRODUCTS • Chemblend Developers • SEF<sup>2</sup> Developer • Autochem® Dental • HS Manual Solutions  
• CBM20 Mixing Station • Imaging Services • Equipment and Accessories